

Precast Gels Plus (Tris-Glycine) compatible with Biorad Miniprotean Electrophoresis Tank

Introduction:

Our Plastic Precast Gels Plus (Tris-Gly) are precast polyacrylamide gels with high-performance and better band resolution.

- > Automated gel casting technology can ensure stability and repeatability.
- Coated plastic plate can effectively reduce the non-specific adsorption of protein and make the protein band sharper and clearer.
- Gel cassettes can be easily opened with the Gel Opener.
- > Precast Gels Plus do not contain SDS. Compatible with denaturing and native electrophoresis.
- Compatible with mainstream mini electrophoresis tanks on the market, such as Bio-Rad,
- Precast Gels Plus are available in different acrylamide concentrations, including gradient 4-15%, 4-20% of separating gels.

Product Details:

- Height of 4% stacking gel: 1.5cm
- Ratio of acrylamide to methylene bis acrylamide: 29:1
- Gel thickness: 1.0mm
- Well format: 10wells/15wells
- Max sample loading volume: 30µL/well for 15wells, 50µL/well for 10wells
- Cassette size: 100×89×4.8 mm
- Gel size: 84×74×1 mm
- Storage: stored at 4°C for 12 months; do not freeze, room temperature: 3 months

POLYACRYLAMIDE PRECAST GELS (TRIS GLYCINE) Suitable for Biorad Electrophoresis Machines.

CAT. #	GRADIENT	Wells	Max volume	Running buffer	Separation Range	Voltage
PSG2001-415	4-15%	10	50µL	Tris-Glycine	20-200kDa	180V
PSG2001-415	4-15%	15	30µL	Tris-Glycine	20-200kDa	180V
PSG2001-420	4-20%	10	50µL	Tris-Glycine	5-200kDa	180V
PSG2001-420	4-20%	15	30µL	Tris-Glycine	5-200kDa	180V



PRE CAST Gel Protocol:

Precast Gels Plus (Tris-Glycine) do not contain SDS and can be used for denaturing and non-denaturing electrophoresis depending on the running buffer.

A.For denaturing electrophoresis:

- 1. Prepare sample: Mix protein sample loading buffer for denaturing electrophoresis with sample in 1:4 ratio (volume) and then heat the samples at 95°C for 5-10 minutes for optimal results.
- 2. Prepare running buffer for denaturing electrophoresis: Use **Tris-Glycine Running Buffer** for **SDS-PAGE** if you have 10× running buffer must be diluted to 1× final concentration before use.

Note: The buffer system of Precast Gels Plus (Tris-Gly) is neutral buffer system.

3. Take the Precast Gels Plus (Tris-Gly) out of the bag and **remove the sealing tape from the top to the bottom of gel cassette**. Assemble the corresponding electrophoresis tank, add the running buffer, and then gently pull the comb out of the cassette.

4. Rinse the wells several times with 1 × running buffer to remove storage buffer before loadingsamples. Load the appropriate volume of your samples in the appropriate wells. Load your protein ladder in the appropriate well.

Note: Do not insert the pipette tip too much into the sample wells to avoid sample leakage caused by deformation of the plastic plates.5. Run the gel: The electrophoresis conditions are usually 180V, 50~60 min.

Note: The voltage is set at 180V constant with an initial current of 75 mA per gel, and a run time of about 60 minutes, and the current will gradually decrease with time.

6. After electrophoresis is complete, retrieve the gel from gel cassette: Inserting the Gel Opener between the two plates and push down gently on the knife handle to separate the plates; repeat on each side of the cassette until the plates are completely separated; then remove the top plate and gently peel away the gel from the plate.



Compatible electrophoresis tanks:

Precast Gels Plus (Tris-Gly) is compatible with most common mini SDS-PAGE tanks, including:

• Bio-Rad Mini-PROTEAN (a simple modification is required, see below)

Precast Gels Plus (Tris-Gly) used in Bio-Rad electrophoresis tank:

- a. Pull out the U-shaped silicone gasket (green), electrode inner core silicone.
- b. Rotate the gasket 180° degrees and install it back oppositely (outside inwards, inside-out). Press down around with your thumbs using even pressure to ensure it seal well to prevent leakage.
- c. Now the system is ready to use.



Pull out the U-shaped silicone gasket

Install the gasket back oppositely

- a. Pull out the U-shaped silicone gasket (green), electrode inner core Silicone (picture A.)
- b. Rotate the gasket 180 degree and Installed it back oppositely (outside inwards, inside-out) picture B.
- c. Press down around with your thumbs using even pressure to ensure to seal well to prevent leakage.
- d. Now the system is ready to use



Distributed by:



Via Rio nell'Elba, 140 – Rome / ITALY Phone +39 06-8818936 / 8800211 – Fax +39 06-8815319 Email <u>info@sichim.com</u> Web-Site <u>www.sichim.com</u>