

ISPKH67 PHK67 Green Fluorescent Cell Linker Kit (in vivo cells, in vitro cells, fixed cells)

Kit components

CAT.#	Description	PKH67 Dye	Diluent C
ISPKH67-100	PKH67 Green Fluorescent Cell Linker Kit Mini	0.1mL	10mL
ISPKH67- 1	PKH67 Green Fluorescent Cell Linker Kit Maxi	1 mL	60mL

Product Description


The PKH67 Fluorescent Cell Linker Kits use proprietary membrane labeling technology to stably incorporate a green fluorescent dye with long aliphatic tails (PKH67) into lipid regions of the cell membrane.^{1,2} The labeling vehicle provided in the kits (Diluent C) is an aqueous solution designed to maintain cell viability, while maximizing dye solubility and staining efficiency during the labeling step.³ Diluent C is iso-osmotic for mammalian cells and contains no detergents or organic solvents, but also lacks physiologic salts and buffers. The appearance of labeled cells may vary from bright and uniform to punctate or patchy, depending on the cell type being labeled and the extent to which membrane internalization occurs after labeling.²⁻⁵ However, PKH67 fluorescence is independent of pH within physiologic ranges and fluorescence intensity per cell is typically unaffected by the pattern of dye localization.

PKH67 fluoresces in the green (Figure 1). Due to its longer aliphatic carbon tails, PKH67 exhibits reduced cell-cell transfer compared with its predecessor, PKH67 is well suited for cytotoxicity assays that use propidium iodide or 7-aminoactinomycin D as viability probes for use in combination with orange-red fluorescent probes such as phycoerythrin, red fluorescent proteins, etc. PKH67 is often used for proliferation monitoring based on dye dilution including estimation of antigen-specific precursor frequencies and identification of quiescent/slowly dividing tumor cells with stem-like properties.

It has also proven useful for monitoring exosome or liposome uptake, cell-cell membrane transfer, phagocytosis, and antigen presentation as well as for **in vivo cell trafficking studies**.

Correlation of in vitro cell membrane retention with in vivo rate of intensity decrease in non-dividing cells predicts an in vivo fluorescence half-life for PKH67 of 10–12 days.

PKH67 is therefore recommended for short-to-medium term in vivo studies requiring a green cell linker dye, as well as for in vitro cytotoxicity, phagocytosis, proliferation, antigen presentation, or other co-culture assays. Due to its extremely stable fluorescence, PKH26 remains the cell linker dye of choice for longer term in vivo studies in which labeled cells are to be followed for periods longer than a few weeks.



Storage and Stability

The PKH67 ethanolic dye solution may be stored at room temperature or refrigerated. To prevent increases in dye concentration due to evaporation, keep the ethanolic dye solution tightly capped except when in immediate use. The dye solution must be protected from bright direct light and examined for crystals prior to use. If crystals are noted in the dye solution, warm slightly in a 37°C water bath, and sonicate or vortex until redissolved. Diluent C may be stored at room temperature or refrigerated. If refrigerated, bring to room temperature before preparing cell and dye suspensions for labeling. Diluent C is provided as a sterile solution. Because it does not contain any preservatives or antibiotics, it should be kept sterile. Do not store dye in Diluent C.

Working solutions of dye in Diluent C should be made immediately prior to use.

Equipment and Reagents

Required but Not Provided for General Membrane Labeling

- Uniform suspension of well-dispersed single cells in tissue culture medium
- Tissue culture medium with serum (complete medium)
- Ca²⁺, Mg²⁺, and serum free medium or buffered salt solution (e.g., Dulbecco's PBS or Hank's BSS)
- Serum, albumin, or other system-compatible protein
- Polypropylene conical bottom centrifuge tubes (4–15 mL)
- Temperature-controlled centrifuge (up to 1,000 g)
- Instrument(s) for fluorescence analysis (fluorescence plate reader, fluorescence or confocal microscope, flow cytometer)
- Laminar flow hood
- Hemocytometer or cell counter
- Slides and coverslip

Procedures

A. General Cell Membrane Labeling.

Labeling occurs by partitioning of the lipophilic dye into cell membranes. Labeling intensity is a function of both dye concentration and cell concentration and is not saturable. Therefore, it is essential that the amount of dye available for incorporation be limited. Over-labeling of cells will result in loss of membrane integrity and reduced cell recovery.

The following labeling procedure can be used for in vitro or ex vivo labeling of stem cells, lymphocytes, monocytes, endothelial cells, neurons or any other cell type where partitioning of dye into lipid regions of the cell membrane is desired. Modified procedures may be required for in vivo labeling and for labeling platelets or phagocytes.

General cell membrane labeling should be performed prior to monoclonal antibody staining. The membrane dyes will remain stable during the monoclonal staining at 4 °C; however, capping of the monoclonal antibodies is highly probable if the general cell membrane labeling is carried out at ambient temperature subsequent to antibody labeling.

The cell and dye concentrations given in the following procedure represent starting concentrations that have been found broadly applicable to a variety of cell types.⁷ Users must determine the optimal dye and cell concentrations for their cell type(s) and experimental purposes by evaluating post-staining cell viability (e.g., propidium iodide exclusion), fluorescence intensity, staining homogeneity, and lack of effect on cell function(s) of interest.

Note 1 : *No azide or metabolic poisons should be present at the time of staining with PKH67.*

Note 2: *Although adherent cells may be labeled while attached to a substrate, more homogeneous staining is obtained using single cell suspensions. Best results will be obtained if adherent or bound cells are dispersed into a single cell suspension using proteolytic enzymes, e.g., trypsin/EDTA, prior to staining.*

The following procedure uses a 2 mL of final staining volume containing final concentrations of 2×10^{-6} M of PKH67 and 1×10^7 cells/mL.

Perform all further steps at ambient temperature (20–25 °C)

1. Place a suspension containing 2×10^7 single cells in a conical bottom polypropylene tube and wash once using medium without serum.

Note: Serum proteins and lipids also bind the dye, reducing the effective concentration available for membrane labeling. Best results are obtained by washing once with serum-free medium or buffer (step 2) prior to resuspension in Diluent C for labeling (step 4).

2. Centrifuge the cells (400 x g) for 5 minutes into a loose pellet.

Note: The PKH67 ethanolic dye solution should not be added directly to the cell pellet. This will result in heterogeneous staining and reduced cell viability.

3. After centrifuging cells, carefully aspirate the supernatant, being careful not to remove any cells but leaving no more than 25 mL of supernatant.

Note: For reproducible results, it is important to minimize the amount of residual medium or buffer present when cells are resuspended in Diluent C. See Note 28 in Ref. 10 for suggested methods.

4. Prepare a 2x Cell Suspension by adding 1 mL of Diluent C to the cell pellet and resuspend with gentle pipetting to insure complete dispersion. Do not vortex and do not let cells stand in Diluent C for long periods of time.

Note: The presence of physiologic salts causes the dye to form micelles and substantially reduces staining efficiency. Therefore, it is important that the cells be suspended in Diluent C at the time dye is added, not in medium or buffered salt solutions.

5. Immediately prior to staining, prepare a 2x dye solution (4×10^{-6} M) in Diluent C by adding 4 μ l of

PKH67 ethanolic dye solution to 1 mL of Diluent C in a polypropylene centrifuge tube and mix well to disperse. **Note 1:** To minimize ethanol effects on cell viability, the volume of dye added in step 5 should result in no more than 1–2% ethanol at the end of step 6.

Note 2: If a final dye concentration $< 2 \times 10^{-6}$ M is desired, the most reproducible results will be obtained by diluting the PKH67 ethanolic dye solution provided in the kit with 100% ethanol to make an intermediate dye stock.

6. Rapidly add the 1 mL of 2x Cell Suspension (step 4) to 1 mL of 2x Dye Solution (step 5) and immediately mix the sample by pipetting. Final concentrations after mixing the indicated volumes will be 1×10^7 cells/mL and 2×10^{-6} M PKH67.

Note 1: Because staining is nearly instantaneous, rapid and homogeneous dispersion of cells in dye solution is essential for bright, uniform and reproducible labeling.

The following measures have been found to aid in optimizing results:

- Do not add ethanolic PKH67 dye directly to the 2x Cell Suspension in Diluent C.
- Mix equal volumes of 2x Cell Suspension (step 4) and 2x Dye Solution (step 5).
- Adjust 2x cell and 2x dye concentrations to avoid staining in very small ($< 100 \mu$ L) or very large (> 5 mL) volumes.
- Use a Pipetman or equivalent for rapid addition of cells and admixing with dye. Serological pipettes are slower and give less uniform staining. Mixing by "racking" or vortexing is also slower and gives less uniform staining.
- Dispense volumes as precisely as possible in order to accurately reproduce both cell and dye concentrations from sample to sample and study to study.

7. Incubate the cell/dye suspension from step 6 for 1–5 minutes with periodic mixing. Because staining is so rapid, longer times provide no advantage.

Note: Expose cells to dye solution and Diluent C for the minimum time needed to achieve the desired staining intensity. Since Diluent C lacks physiologic salts, longer exposure may cause reduced viability in some cell types. If such effects are suspected, include a diluent-only control and a mock-stained control using ethanol rather than dye.

8. Stop the staining by adding an equal volume (2 mL) of serum or other suitable protein solution (e.g., 1% BSA) and incubate for 1 minute to allow binding of excess dye.

Note 1: Serum (or an equivalent protein concentration) is preferred as the stop solution. Increase volume to 10 mL if complete medium is used instead of serum.

Note 2: Do not stop by adding Diluent C or centrifuge the cells in Diluent C before stopping the staining reaction.

Note 3: Do not use serum-free medium or buffered salt solutions, which cause formation of cell-associated dye aggregates. Dye aggregates act as slow-release reservoirs of unbound dye that are not efficiently removed by washing and can transfer to unlabeled cells present in an assay.

9. Centrifuge the cells at 400xg for 10 minutes at 20-25 °C and carefully remove the supernatant, being sure not to remove cells. Resuspend cell pellet in 10 mL of complete medium, transfer to a fresh sterile conical poly-propylene tube, centrifuge at 400xg for 5 minutes at 20-25 °C, and wash the cell pellet 2 more times with 10 mL of complete medium to ensure removal of unbound dye.

Note 1: Transfer to a fresh tube increases washing efficiency by minimizing carryover of residual dye bound to tube walls.


Note 2: Do not use Diluent C for washing steps.

10. After the final wash, resuspend the cell pellet in 10 mL of complete medium for assessment of cell recovery, cell viability and fluorescence intensity (Figure 2). Centrifuge and resuspend to desired final concentration of viable cells.

Note 1: Stained cells may be fixed with 1-2% neutral buffered formaldehyde and intensities are stable for at least 3 weeks if samples are protected from light.

Note 2: Staining is typically at least 100–1,000 times brighter than background autofluorescence. Intensity distributions should be symmetrical and as homogeneous as possible, although staining

Figure 2. Staining Optimization for PKH67




PKH67 staining concentration may be optimized using a method similar to that described for PKH26.35 MC-38 TIL cells were stained with the indicated concentrations of PKH26 dye at a final cell concentration of 1×10^7 cells/ml. Viability (\blacktriangle) was determined by dye exclusion and mean fluorescence intensity (\bullet) was determined by flow cytometry. Anti-tumor TIL specificity and potency in vitro and in vivo was confirmed to be unaltered by labeling with 20 μ M PKH26.35

B. Histology Fixed cells and tissues

PKH26 labeled adipocytes have been successfully identified in tissues subjected to standard paraffin embedding and sectioning, but such methods risk loss of intensity because they use clearing agents that may partially extract membrane lipids and lipophilic dyes. Histologic studies of tissues containing cells labeled with lipophilic membrane dyes have typically been carried out on serial frozen sections or sections prepared after fixation in neutral buffered formalin prepared from 4% paraformaldehyde. These methods avoid quenching of fluorescence by absorbing dyes found in histologic counterstains. For studies where imaging is to be done on a single section, fluorescence microscopy should precede counterstaining.

Fluorescence imaging of frozen sections:

1. Remove coverslips by soaking slides in acetone for 24–48 hours.
2. Rinse slides in distilled water to remove acetone.
3. Counterstain sections using stain of choice. Satisfactory results have been obtained using Mayer's or Harris
4. Mount slides using AS/AP permanent aqueous mounting medium.



Internalization kinetics of PKH67 in L929 cells.

Dye	PKH67
Emission	GREEN (max. 502nm)
Useful laser lines	457nm, 488nm
Spectrally compatible with ²	Hoechst 33342, PE, Cy 3, TR, PI, some red FPs, PE-Cy5, PerCP, 7-AAD, APC, TO-PRO-3, DRAQ5, PE-Cy7, APC-Cy7
NOT suitable for use with ²	FITC, CFSE, green/yellow FPs

**For Research use only
IMMUNOLOGICAL SCIENCES**

Web-site: www.immunologicalsciences.com - E-Mail: info@immunologicalsciences.com