


ISB0026 7-AAD Viability Staining Solution

Size:	400 tests
Reactivity:	Broad
Application:	FC (QC tested)
Description:	The 7-AAD (7-aminoactinomycin D) Staining Solution is intended to evaluate cell viability by flow cytometry.
Excitation laser:	blue (488 nm)
Formulation:	7-aminoactinomycin D (7-AAD) diluted in PBS (0.05 mg/ml) containing 15 mM sodium azide.
Storage:	Store in the dark at 2-8°C. Do not freeze. Avoid prolonged exposure to light. Do not use after expiration date stamped on vial label.



Flow cytometry staining of apoptotic JURKAT cells with 7-AAD and annexin V FITC.